




SB-SHIFIX50 COVERSGLIPS

SKU: SB-SHIFIX50



DESCRIPTION

SHI-FIX™ COVERSGLIPS

Do you wish immunofluorescence of suspension cells was as easy as with adherent cells? Now it is! Shi-fix™ (manufactured by Shikhar Biotech) cover-slips allow suspension cells to be layered on them or directly grown as a monolayer. Easy, hassle-free protocol and no spin columns needed for fixing suspension cells to slides. Just add your cells to the coverslips, wait for 30 mins, wash unbound cells with PBS and proceed to immunostaining (or continue the culture to obtain desired cell densities).

These cover-slips can also be used for staining suspension cells for microscopy, or for fixing suspension cells for confocal microscopy. What's more, you can also stimulate the cells while they are attached on the cover-slip, should you need to do so. You can work on non-adherent cells as if they were adherent!



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

Suspension cell immunofluorescence with Shi-fix™ coverslips

Size: ø 12mm x 0.17mm(#1.5)

SPECIFICATIONS

Unit Size

Shipping Instructions

Pack of 50

Ambient

GALLERY IMAGES

