


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-AP3B1 / HPS ANTIBODY





SPECIFICATIONS

Unit Size 100 μg

Synonym / Alias

 $AP3B1|PE|HPS|HPS2|ADTB3|ADTB3A|adaptor-related\ protein\ complex\ 3,\ beta\ 1\ subunit|beta3A-adaptin|AP-3B-adaptin|AP-3B-adaptin|AP-3B-adaptin|AP-3B-adaptin|AP-3B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|$

complex beta-3A subunit

Names Usage

Flow Cytometry: Flow cytometric analysis of HepG2 cells. Recommended concentration:

Summary

10ug/ml.

Accession ID

NP_003655

Blocking

EBP06382

Peptide

Immunogen

 $Peptide \ with \ sequence \ SSNSFPYNEQSGG-C, \ from \ the \ N \ Terminus \ of \ the \ protein \ sequence \ according \ to$

NP_003655.

Peptide Sequence

SSNSFPYNEQSGG-C





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Shipping Refrigerated Instructions

Predicted

Human **Species**

Reactive **Species**

Human

Human **Gene ID** 8546

Product Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite_medium.png$

ELISA

Detection Antibody detection limit dilution 1:64000.

Limit

 $\begin{array}{c} \textbf{Application} \\ \textbf{Pep-ELISA, FC} \end{array}$

Type

DOCUMENTS

• Data Sheet

GALLERY IMAGES

